WINKELSTABILE PROXIMALE TIBIAPLATTE


Fixed-angle proximal tibia plate

- T-Design

GROSSFRAGMENT

- Big fragment

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Fractures of the proximal tibia are common injuries of the lower limb. Plates are the implant of first choice for internal fracture stabilisation particularly when the joint surfaces are involved. The aim of operative management is secure fragment fixation by the implant. After operative restoration of the anatomy, the stability of the treatment should enable early functional follow-up treatment. The principle of locking plate internal fixation offers the familiar advantages of an internal fixator in implant fixation in the joint-bearing part of the bone. In elderly patients in whom bone quality is often osteoporotic, the danger of implant loosening is reduced. The proximal tibial plate is fixed laterally. It is available in different lengths and is anatomically shaped. This shortens the operating time considerably as time-consuming shaping of the implant on the operating table can be omitted. Ideally the plate can also be used as a template for restoring the shape of the bone. The screws in the joint block of the bone are inserted to provide angular stability. In order to achieve optimal force distribution in the bone and prevent secondary sinking of the impression fragments, which are usually situated laterally, three screws can be inserted in the proximal transverse limb of the plate. In the plate shaft in the distal direction, it is possible for you to choose between conventional and locking screw placement in the plate because of the special design of the screw holes (combi-hole). With angularly stable or locking insertion, the screws diverge alternately. This way notching tensions in the bone and the risk of secondary implant migration are diminished.

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- Fractures of the proximal Tibia
- Multi-fragment fractures, metaphyseal
- Joint fractures, inclusive C3 fractures

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

Place the patient supine with a thigh tourniquet and the leg freely mobile. Ensure that the operation site can be shown intraoperatively on the image intensifier both a. p. and laterally. Select a lateral approach; dissect the muscles gently subfascially. The dissection must extend proximal to the fibulo-tibial joint, if necessary with exposure of the peroneal nerve. Expose the site distally to the shaft according to the planned plate length; with the „pushed-in„ technique, a smaller approach is needed. The first stage of internal fixation consists of anatomical reduction and retention of the joint fragments.

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
Select a plate of suitable length and place it strictly laterally on the tibial head until an optimal fit with the proximal surface of the bone is achieved. Please note that there is a right and left version. And the tip of the transverse cranial limb of the plate must point anteriorly. The plate must not cross the level of the joint and must not compromise the fibulo-tibial joint. In the metaphyseal section of the plate, alterations in shape are usually not necessary. With long plates, the individual varus/ valgus and torsion of the tibia must be taken into account and the distal limb of the plate must be shaped appropriately if necessary. Then fix the plate to the bone.

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
Check the result of reduction clinically and radiologically (a.p. and lateral).

In the case of metaphyseal comminuted zones and fractures extending to the shaft, rotation with 90° flexion of the hip and knee and the axis should be checked with the cable method. Fill the proximal angle-stable holes in the plate first. Screw the drill sleeves into the holes, drill with the 3.2 mm drill and after measuring the length with the gauge and removing the drill sleeves insert the 6.0 mm angle-stable screws tightly into the plate bearing. The screws should not cross the medial cortex. The screws fixed in the transverse proximal limb of the plate converge. The subsequent distal screws point cranially and medially. This optimal screw configuration of the joint block enables unilateral stabilisation even with bicondylar fractures.

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
After a further clinical and radiological check, the joint block is reduced definitively to the shaft and retained. With metaphyseal fractures and fractures running distally to the shaft, the rotation and axis are again checked critically. Now fill the remaining screw holes. At least three screws should be inserted bicortically distal to the fracture. With angle-stable fixation the drill sleeve is first screwed into the plate hole, the 3.2 mm drill is used and after measuring the length, the self-tapping threaded-head cortical screw is inserted, alternatively 4.5 mm cortical screws can be used. At the conclusion of the internal fixation, perform a clinical check and radiological check in 2 planes. Deep drains are placed and the wound is closed in layers. Early functional follow-up treatment, limited by the fracture pathology, is usually possible.

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Cancellous screws, D 6,0 mm (conical head thread, self-tapping)

Bestell-Nr. Stahl Code N° steel	Bestell-Nr. Titan Code N° titanium	Länge Length
2.163.65	3.163.65	65 mm
2.163.70	3.163.70	70 mm
2.163.75	3.163.75	75 mm
2.163.80	3.163.80	80 mm
2.163.85	3.163.85	85 mm
2.163.90	3.163.90	90 mm

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

No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying, recording or otherwise without the prior written permission of the publisher.

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© Copyright Königsee Implantate GmbH
OSTEOSYNTHESE GmbH
KÖNIGSEE IMPLANTATE UND INSTRUMENTE ZUR OSTEOSYNTHESE GmbH