Winkelstabile Distale Femurplatte

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Angle-stable femoral plate

Large fragment

Stahl / Titan

Stainless steel / Titanium

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Distal femoral fractures often present the surgeon with a great challenge. Whereas younger patients often sustain complex fractures involving the joint with major soft tissue injury due to high-speed trauma, in geriatric traumatology the problem is one of osteoporotic fractures with reduced bone substance. Furthermore, periprosthetic femoral fractures will become increasingly important in future. Both as regards the fracture and also because of the accompanying circumstances, implants are necessary that ensure the best possible stability and that can also be inserted minimally invasively, thus sparing the soft tissue. The distal femoral plate from Königsee meets these requirements.

The design was tested in an anatomical study in cadaver bone and using conventional radiology. The plate is anatomically shaped. On the other hand, it is not applied directly to the bone so that the periosteal circulation can be preserved. Fixed-angle interlocking is possible at all of the holes in the plate. By using the fluoroscopically visible guide frame, the screws can be placed in the holes in the shaft region without complex access. In the case of periprosthetic fractures with a prosthesis stem, stabilisation can be achieved by means of short screws inserted through the guide frame.

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
The distal femoral plate is employed in the management of fractures of the distal femoral shaft and fractures close to or involving the joint:

- Condyle fractures with joint involvement
- Supracondylar fractures
- Distal femoral shaft fractures
- Periprosthetic femoral fractures

The plate is particularly suitable for management of 33A fractures and 33C fractures.

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
- female, 38 years, road traffic accident
- AO: 33A3-fracture

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The postoperative treatment should be determined by the surgeon.

Primary full loading will usually not be achievable. However, stable movements should always be achieved with treatment on the motor rail. An increase to full loading is to set individually.

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
Use of planning templates assists preoperative planning. These allow the plate length to be determined before the start of the surgery. The operation is performed with the patient supine on an X-ray permeable table. If the fracture is located very far distally, it is sufficient to lower the healthy leg a little. More proximal fractures require positioning the healthy leg in a leg holder so that X-rays can be taken in the lateral plane also. It should be operated in order to minimise the muscle pull on the condyles, which facilitates reduction. After positioning, frontal and lateral X-rays should be taken again before the start of the operation. The entire region of the plate being implanted must be imaged.

A tourniquet can be used only in the case of smaller accesses to manage entirely distal fractures.

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
To treat distal femoral fractures without joint involvement, the lateral access should be selected.

The skin incision is longitudinal and the iliotibial tract is split. The vastus lateralis muscle is held forward by a Hohmann retractor. This means that the outside of the lateral femoral condyle is exposed. Intraarticular fractures with joint

involvement require anterolateral access with arthrotomy. This allows reduction and retention under vision.

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Apart from the very distal fractures, which allow use of the 4-hole plate, the guide frame should be used for soft tissue sparing implantation.

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
The distal part of the guide frame fits the plate design and sits directly on the plate.

Plates and guide frames are identified with "R" and "L" for the corresponding side.

A long threaded drill bush for 3.2 mm drill is first inserted through the most distal hole of the projecting part of the guide frame and fixed in the plate.

To tighten it (caution: do not damage the guide frame), the turn pin can also be used. A retaining cap is then screwed on to ensure stability. A drill bush for Kirschner wires is then inserted through the guide attachment into the most distal hole of the part of the guide frame in contact with the plate. The plate is thus fixed to the guide frame with adequate stability.

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

In the management of fractures involving the joint, reconstruction of the articular fragment is necessary first. This can be done indirectly under fluoroscopic control or directly by opening the joint. Depending on the fracture situation, retention can be temporary with Kirschner wires or definitive with screws. The screws should be positioned so that they do not interfere with subsequent implantation of the plate. In the metaphyseal region, reduction is facilitated by lowering the lower leg. The fracture can usually be aligned indirectly by pushing the plate proximally. In the fracture region, reduction is possible in the form of the Rep-ret-system. The threaded wire is inserted through the drill bush. Compression can be obtained with the wing nut through the retaining cap.

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
Following dissection of the lateral femoral condyle, the plate with attached guide frame can be inserted. It is useful to prepare the way between the vastus lateralis and periosteum with the tissue dissector. The plate with guide frame can then be pushed in. It should first be guided further proximally and then back and distally until it sits correctly on the lateral condyle.

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
A stab incision is now made at the level of the proximal end of the plate. The plate lengths are clearly marked on the guide frame. Depending on the length of the selected plate, a threaded tissue protection sleeve containing a trocar is inserted. This is why the guide frame holes for the most proximal plate holes of all plate lengths are threaded.

After removing the trocar, a drill bush for 1.6 mm Kirschner wires is screwed in and then stabilised with a retaining cap. A Kirschner wire is then screwed through both cortices. The proximal plate end must sit on the centre of the femoral shaft. The position should be checked by a lateral X-ray.

The system is thus closed and the plate is fixed to the femur with adequate stability.

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For fixation of the plate it is important that it is not in a tilted position on the femur. Because of the intrinsic weight of the guide frame, the plate-guide-frame construction has a tendency to tilt away from the femoral shaft. This can be prevented by carefully pressing the distal end of the plate under vision. Through the distally inserted drill bush for Kirschner wires, a 1.6 mm wire is then drilled. This establishes the distal plate position.

Ensure by X-ray that this wire is at least 3 mm from the joint so that the screw that will be inserted later does not injure the cartilage or enter the joint. Optionally, two more Kirschner wires can also be inserted through the small holes at the distal end of the guide frame.

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The screws can now be inserted according to the preoperative planning. The screws should first be placed in the condylar region according to the part of the guide frame in contact with the plate. Selftapping fixed-angle 5.5 mm cortical screws are available.

The tissue protection sleeve is inserted and the drill bush for 4.0 mm spiral drill is screwed into it.

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A preliminary hole is drilled with the 4.0 mm spiral drill. The length is measured by means of the advanced measuring aid.

Alternatively, the length can also be determined with a gauge.

To do so, the drill bush is removed and the gauge is put on directly on the plate hole.

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

If non-fixed-angle management is planned distally, the drill bush for 3.2 mm spiral drill must be used as well as the tissue protection sleeve. The 3.2 mm drill is used for drilling.

For fixed-angle screws in general: measurement reading = screw length. 2 mm must be added for non-fixed-angle screws. Finally, the screw is inserted with the hexagonal screwdriver. Screws are inserted the same way in all the other distal plate holes.

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For filling the shaft holes, following a stab incision the unthreaded tissue protection sleeve is advanced to the plate together with the trocar through the corresponding hole in the guide frame. After removing the trocar, the drill bush is screwed in. A hole is predrilled with the 3.2 mm drill and the attached measurement aid is used to determine screw length.

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

Alternatively, the length can also be measured with the gauge after removing the drill bush. The distance sleeve remains on the gauge. Drilling is bicortical (interior implants are exceptions). The cortical screw is screwed in through the tissue protection sleeve with the hexagonal screwdriver. Finally, the filled plate hole is identified by snapping in a blind plug.

All of the other proximal shaft holes are filled in the same way. If non-fixed-angle cortical screws are used, the screws should be 2 mm longer. When fixed-angle screws are used, at least 4 screws should be placed in the shaft above the fracture. If necessary, any plate holes not occupied by screws can be filled with a screw head to prevent tissue migration into the plate.

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
After all the necessary shaft screws have been inserted, the guide frame must be removed before the last screw is placed on the condyle.

The 4.0 mm drill can be used to drill through the drill guide and the screw can then be inserted.

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
After placing a deep drain, the wound is closed in layers with subcutaneous drain and skin suture.

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To remove the implant after the fracture has healed, the guide frame should be used with longer plates in order to minimise soft tissue trauma.

The proximal screws can be removed through the inserted sleeves and the distal ones can be unscrewed directly.

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Instruments in extension-set, steel and titanium

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1

Gewebelöser

Instrument for taking off the tissue

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10.176.22

1

Bohrbuchse Drill guide for K-wire D 1.6 mm, lenght 40 mm, ID 1.8 mm, AD 10.0 mm

für K-Draht, D 1,6 mm, Länge 40 mm, ID 1,8 mm, AD 10,0 mm

10.176.20

1


Bohrbuchse Drill guide for K-wire D 1.6 mm, lenght 171 mm, ID 1.8 mm, AD 8.0 mm

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Locating sieve (for instruments extension-set distal femur)

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Hexagonal screw-driver, with handle, for screws D 4,5 mm – 7,0 mm, length 260 mm

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Sleeve for femur guide frame, D 10.0 mm, long, with thread

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Sleeve for femur guide frame, D 10.0 mm, long, without thread

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Trocar, to transfix soft tissues, D 8.0 mm

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D 3,2 mm, Länge 305 mm

Twist drill (for quick coupler) D 3.2 mm, length 305 mm

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D 4,0 mm, Länge 305 mm

Twist drill (for quick coupler) D 4.0 mm, length 305 mm

10.261.14F

1

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Measure instrument, for drill D 3.2 mm

2.952.032

2

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Measure instrument, for drill D 4.0 mm

2.952.040

2

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Screw gauge, with clasp for fixed angle screws, MB 110 mm

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Kirschner wire with trocar point and round end, D 1.6 mm, length 280 mm

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Cortical screws, diameter 4,5mm (self-tapping)

2.152.14	3.152.14	14 mm	6
2.152.18	3.152.18	18 mm	6
2.152.22	3.152.22	22 mm	6
2.152.26	3.152.26	26 mm	6
2.152.30	3.152.30	30 mm	6
2.152.34	3.152.34	34 mm	6
2.152.38	3.152.38	38 mm	6
2.152.42	3.152.42	42 mm	6
2.152.46	3.152.46	46 mm	6
2.152.50	3.152.50	50 mm	3
2.152.54	3.152.54	54 mm	3

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Cortical screws, diameter 4,5mm, fixed angle, conical head thread (self-tapping)

2.545.14	3.545.14	14 mm	6
2.545.18	3.545.18	18 mm	6
2.545.22	3.545.22	22 mm	6
2.545.26	3.545.26	26 mm	6
2.545.30	3.545.30	30 mm	6
2.545.34	3.545.34	34 mm	6
2.545.38	3.545.38	38 mm	6
2.545.42	3.545.42	42 mm	6
2.545.46	3.545.46	46 mm	6
2.545.50	3.545.50	50 mm	3
2.545.54	3.545.54	54 mm	3

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Cortical screws, diameter 5,5mm, fixed angle, conical head thread (self-tapping)


2.549.40	3.549.40	40 mm	2
2.549.45	3.549.45	45 mm	2
2.549.50	3.549.50	50 mm	2
2.549.55	3.549.55	55 mm	4
2.549.60	3.549.60	60 mm	4
2.549.65	3.549.65	65 mm	4
2.549.70	3.549.70	70 mm	4
2.549.75	3.549.75	75 mm	4
2.549.80	3.549.80	80 mm	4
2.549.85	3.549.85	85 mm	4
2.549.90	3.549.90	90 mm	2

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Screw head, conical thread (hole plug)

8.946.05	9.946.05		5
----------	----------	--	---

OP-Set Stahl OP-Set, steel

Bestell-Nr. 19.600.00
Code No

OP-Set Titan OP-Set, titanium

Bestell-Nr. 19.602.00
Code No

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Instruments in set, steel and titanium

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Locating sieve (for instruments and implants)

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Hexagonal screw-driver, with handle, for screws D 4.5 mm – 7.0 mm

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for angle stable screwing D 4.5 mm, for twist drill D 3.2 mm

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for angle stable screwing D 5.5 mm, for twist drill D 4.0 mm

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D 3,2 mm, Länge 195 mm

Twist drill (for quick coupler) diameter 3.2 mm, length 195 mm

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D 4,0 mm, Länge 220 mm

Twist drill (for quick coupler) diameter 4.0 mm, length 220 mm

2.904.23 1

Schraubenhaltepinzette

Screw-holding forceps

2.954.01 1

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Screw gauge, with clasp for fixed angle screws, MB 90 mm

2.953.90 1

Kirschnerdraht Trokar, Ende rund, D 1,6 mm, Länge 150 mm

Kirschner wire with trocar point and round end, D 1.6 mm, length 150 mm

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
Angle stable distal femur plate, right, for screws D 4.5 – 5.5 mm

Bestell-Nr. Stahl Code No steel	Bestell-Nr. Titan Code No titanium	Länge Length	Lochanzahl No of hole	Anzahl im Set Quantity in set
4.9462.06	5.9462.06	150 mm	7/6	1
4.9482.08	5.9482.08	182 mm	7/8	1
4.9482.12	5.9482.12	246 mm	7/12	1
4.9482.16	5.9482.16	309 mm	7/16	1

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Angle stable distal femur plate, left, for screws D 4.5 – 6.0 mm

Bestell-Nr. Stahl Code No steel	Bestell-Nr. Titan Code No titanium	Länge Length	Lochanzahl No of hole	Anzahl im Set Quantity in set
4.9472.06	5.9472.06	150 mm	7/6	1
4.9492.08	5.9492.08	182 mm	7/8	1
4.9492.12	5.9492.12	246 mm	7/12	1
4.9492.16	5.9492.16	309 mm	7/16	1

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
*No part of this publication may be reproduced,
stored in a retrieval system or transmitted in any
form or by any means electronic,
mechanical, photocopying, recording or other-
wise without the prior written permission of the
publisher.*

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© Copyright Königsee Implantate GmbH
KÖNIGSEE IMPLANTATE UND INSTRUMENTE ZUR OSTEOSYNTHESE GmbH